

New Instructions for Patient CD

1. Insert the disc into the CD drive.
2. The CD will automatically open to the MRI Group logo.
3. Click on the **View Images** icon in the upper right hand toolbar.
4. This will launch you into the image viewer.
5. Double click on your listed exam or click on exam and touch **View** at the bottom of the page.
6. This will launch your images.
7. The default is set to left mouse to scroll through the images and right mouse to change the contrast.
8. There is also a **Help** button in the right side of the toolbar for additional usages of this image viewer.
9. To exit select **File** in the top left toolbar and then select **Exit**.

Old Instructions for Patient CD

1. Insert disc into CD drive
2. Click "ok" on Syngo Fastview screen
3. Highlight patient's name and "drag" into black background
4. "X" out of the browser

You should now have an image on the left side of the screen and two boxes on the right side of the screen with options for viewing.

To Page Through The Images

1. Click on the "scroll" tab on the bottom box
2. By selecting "patient + and patient -" options on top row, you will be able to page through the images.

To View Images in Different Formats

1. Select the "view" tab on the top box
2. Select format you would like to view images in, for example 1 square, 4 squares or 9 squares

To Change Contrast in Images

1. Select "image" tab on the top box
2. Select the first icon in the box
3. Place your cursor over the image you would like to view
4. Hold the mouse down and move the cursor back and forth. This will change the image contrast.

To Exit

1. Select "file" from the top of the tool bar
2. Exit